

Archtek

Bleaching and Utility Dental Trays



■ **Archtek's Hinged Arch** bleaching tray in 1.0 mm thickness requires no boiling to form (#133)



■ **Archtek's Hot Pot** generates warm water in 2 minutes, always to the exact temperature needed (#695)



■ Only **Archtek** offers four styles of “Boil & Bite” Dental Trays

- **Full-Arch** standard tray in 1.5 mm thickness (#112)
- **Two-Tray System** tray in 1.0 mm thickness (#117)
- **Ultra-Thin Partial Arch** tray in 1.0 mm thickness (#124)
- **Ultra-Thin Full Arch** tray in 1.0 mm thickness (#125)

■ The “Boil & Bite” forming allows for immediate chair-side fitting

■ The front tab makes forming easy and is removed once successfully fitted

■ All trays can be trimmed to the gingival margins once the tray is formed

■ All styles are available in bulk or in convenient kit packs with re-sealable storage bags

Archtek
INNOVATION IN DENTISTRY

12105 West Cedar Dr. | Lakewood, CO 80228 USA

303/763-8916 | 800/763-8129 TOLL-FREE | 800/763-1915 FAX

www.archtekinc.com WEB | sales@archtekinc.com E-MAIL